

Serial No. 09/915,411

PAGE 3
CASE 51882AUSM1SPECIFICATION AMENDMENTS**EXAMPLE 1***In Vitro* Assay*In vitro* binding assay for non-peptide CCR1 receptor antagonists

This assay demonstrates the affinities of the non-peptide CCR1 receptor antagonists of the invention, preferably a non-peptide CCR1 receptor antagonist of formula (I), for binding to the rat CCR1 receptor.

Reagents and Solutions:Chemokines: MIP-1 α and RANTES (Peprotech Inc.)

Cells: Rat peripheral blood mononuclear cells (PBMC) were isolated from whole blood from Lewis rats by Accu-Paque™ ACCU-PAQUE™ (4.2% w/v dextran and 11.6% w/v sodium diatrizoate) (Accurate Chemical & Scientific Corp.) density centrifugation.

Ligand: ¹²⁵I-MIP-1 α and ¹²⁵I-RANTES from New England Nuclear (specific activity is 2200 Ci/mmol, 25 μ Ci/vial) was reconstituted in 1 mL H₂O.

Assay buffer: 130 mM NaCl, 5 mM KCl, 1 mM MnCl₂, 50mM Tris, 30 μ g/ml bacitracin, 0.1% BSA, pH 7.4.

Wash buffer: Phosphate buffer solution (PBS)

Compounds of the Invention: The stock solution of the compounds was 1 mM in 100% DMSO.

The highest concentration in the assay was 10 μ M and may vary depending on the potency of the compounds. Serial 1:3 dilutions from the highest concentration were made with assay buffer. Six concentrations of each compound were typically screened to generate a dose curve from which the K_i value was determined.

Assay procedure:

Assays were performed in 96-well v-bottom microtiter plates in a total volume of 100 μ L.

Rat PBMC were washed once in PBS and resuspended in the assay buffer to about 0.2 to 1.0 x 10⁶ cells/mL. Cells were incubated with ¹²⁵I-MIP-1 α or ¹²⁵I-RANTES in the presence or absence of varying concentrations of unlabeled MIP-1 α , RANTES, or compound at 4°C for 30 minutes.

The reactions were terminated by removing aliquots from the cell suspension and separating cells from buffer by centrifugation through a silicon/paraffin oil mixture as described in Hesselgesser *et al.*, (1998), *supra*.

n-0500

Serial No. 09/915,411

PAGE 4
CASE 51882AUSM1

The nonspecific binding was determined in the presence of 100 nM or 1 μ M of unlabeled MIP-1 α or RANTES. The concentrations of compounds in the assay were typically from 10 μ M to 30 nM in 1:3 dilutions and the concentrations for more potent compounds were lower depending on the potency.

Calculations:

The dose curves of each compound with 6 concentration points were generated and the binding data were curve fitted with the computer program IGOR (Wavemetrics) to determine the affinity and number of sites.

The non-peptide CCR1 receptor antagonists of the invention, when tested in this assay, demonstrated their affinity to bind to the rat CCR1 receptor.

n-0500